

SOC Medium

Catalogue No. AS-1428

Date of issue: 27/09/2024 | Revision: 1.0

Product Information

Product name	SOC Medium (dehydrated base)
Catalogue number	AS-1428
Other names	S.O.C. Medium; SOC Broth; Super Optimal broth with Catabolite repression
Product type	Dehydrated molecular biology culture medium (glucose-free base)
Usage	28.09 g in 980 mL, plus 20 mL of sterile 1 M glucose after autoclaving
Final pH	7.0 ± 0.2 at 25 °C
HS code	3821.00

Intended Use

SOC Medium is used for the outgrowth of *Escherichia coli* after chemical transformation or electroporation. The glucose-free base can also be used as SOB for growing cells before making them competent. For laboratory microbiological testing, quality control and research use only. Not intended for direct use in food, human or veterinary therapeutic applications, diagnosis, or consumption.

Summary and History

Hanahan (1983) developed SOB and SOC while optimising plasmid transformation of *E. coli*. Outgrowth in SOC gives more transformants than outgrowth in LB, and SOC is now the standard recovery medium in most cloning protocols.

Mode of Action

Tryptone and yeast extract supply amino acids, vitamins and nucleotides for rapid recovery. Glucose provides a readily used energy source, and magnesium supports membrane and ribosome function. Low salt limits osmotic stress while cells repair their membranes and express the plasmid resistance gene.

Typical Formula (per litre of final SOC)

Ingredient	g/L	Concentration	Function
Tryptone (pancreatic digest of casein)	20.0	2%	Amino acids and peptides
Yeast extract	5.0	0.5%	Vitamins, nucleotides, growth factors
Sodium chloride	0.5	8.6 mM	Low ionic strength
Potassium chloride	0.186	2.5 mM	Potassium source
Magnesium sulfate (anhydrous)	2.4	20 mM	Magnesium source
Dehydrated base total	28.09		Suspend in 980 mL
Glucose (added after autoclaving, not in powder)	3.6	20 mM	Energy source

Hanahan's original SOC supplies 20 mM Mg²⁺ as 10 mM MgCl₂ plus 10 mM MgSO₄; AS-1428 supplies it as MgSO₄. Formula may be adjusted to meet performance requirements.

Preparation

- Dissolve 28.09 g of powder in 980 mL of purified water.

- Autoclave at 121 °C for 15 minutes.
- Cool to below 50 °C.
- Aseptically add 20 mL of sterile 1 M glucose (18% w/v), sterilised separately by 0.2 µm filtration. Do not autoclave glucose together with the base.
- Mix and dispense aseptically. Final volume: 1 L.

Typical Transformation Outgrowth

- After heat shock or electroporation, add 250 µL to 1 mL of room-temperature SOC to the cells.
- Incubate at 37 °C for 1 hour with shaking at about 225 rpm.
- Spread on selective agar plates and incubate overnight at 37 °C.

Quality Control

Organism / test	Inoculum (CFU)	Incubation	Expected result
<i>Escherichia coli</i> K-12 ATCC 10798	10–100	37 °C, 18–24 h	Good growth
Uninoculated medium (sterility)	—	37 °C, 48 h	No growth

Limitations

- SOC is not selective. Plate onto media containing the appropriate antibiotic after outgrowth.
- The medium is formulated for *E. coli*; it is not intended for yeast transformation recovery.
- Transformation efficiency depends mainly on the competent cells, DNA and protocol.

Storage and Shelf Life

Dehydrated base	Store at 15–25 °C in a dry place with the container tightly closed. Use before the expiry date on the label.
Prepared SOC	Store at 2–8 °C, or in single-use aliquots at –20 °C. Discard if cloudy. Validate the storage period in-house.
Signs of deterioration	Caking, discolouration or moisture absorption in the powder.

References

- Hanahan D. Studies on transformation of *Escherichia coli* with plasmids. *J Mol Biol.* 1983;166(4):557–580.
- Green MR, Sambrook J. *Molecular Cloning: A Laboratory Manual*. 4th ed. Cold Spring Harbor Laboratory Press; 2012.
- Inoue H, Nojima H, Okayama H. High efficiency transformation of *Escherichia coli* with plasmids. *Gene.* 1990;96(1):23–28.

The information in this document is based on current technical knowledge and is provided in good faith for laboratory and technical evaluation purposes. It does not constitute a warranty of fitness for a particular purpose. Users are responsible for validating the product within their own methods and for complying with applicable regulations. AuSaMicS Pty Ltd accepts no liability for loss or damage arising from use of this information. © AuSaMicS Pty Ltd. All rights reserved.