

Ceftriaxone Disodium Salt Hemi Heptahydrate

Catalog No. AS-2007 | CAS 104376-79-6

⚠ For Laboratory Use Only — Not for Human or Veterinary Use

1. Product Identification

Product Name	Ceftriaxone Disodium Salt Hemi Heptahydrate
Synonyms	Disodium Ceftriaxone Salt; Ceftriaxone Na ₂ ; Rocephin® (pharmaceutical grade)
Catalog Number	AS-2007
CAS Number	104376-79-6
Molecular Formula	C ₁₈ H ₁₆ N ₈ Na ₂ O ₇ S ₃ · 3.5H ₂ O
Molecular Weight	661.60 g/mol
Class	Third-Generation Cephalosporin Antibiotic
Grade	Research Grade
Supplier	AuSaMicS Pty Ltd 31 Longview CT, Thomastown VIC 3074 Australia
ABN	56 676 640 467

2. Physical & Chemical Properties

Appearance	White to slightly yellowish crystalline powder
Odour	Odourless
Solubility	Freely soluble in water; slightly soluble in methanol; practically insoluble in ethanol and acetone
pH (1% aq. sol.)	6.0 – 8.0
Melting Point	Decomposes above 155 °C (no defined melting point)
Storage Temp.	2 – 8 °C, protected from light and moisture
Hygroscopicity	Hygroscopic — store in tightly sealed container
Specific Rotation	[α] _D ²⁰ = –150° to –170° (c = 1, H ₂ O)

3. Biological Activity & Mechanism of Action

Ceftriaxone is a potent third-generation cephalosporin antibiotic with broad-spectrum bactericidal activity against both Gram-positive and Gram-negative organisms. It exerts its antibacterial effect by irreversibly inhibiting the

transpeptidase enzymes (penicillin-binding proteins, PBPs) involved in the final cross-linking step of peptidoglycan biosynthesis. This disrupts the structural integrity of the bacterial cell wall, leading to osmotic lysis and cell death.

Target	Penicillin-Binding Proteins (PBPs 1a, 1b, 2, 3)
Action	Inhibition of peptidoglycan cross-linking → cell wall lysis
Gram-positive spectrum	Staphylococcus aureus (MSSA), Streptococcus spp., Streptococcus pneumoniae
Gram-negative spectrum	Escherichia coli, Klebsiella spp., Haemophilus influenzae, Neisseria spp., Proteus mirabilis, Salmonella spp.
MIC Range	0.001 – 0.5 µg/mL (organism-dependent)
Resistance Mechanism	β-lactamase hydrolysis; PBP alteration; efflux pumps

4. Preparation & Reconstitution Guidelines

Recommended Solvent	Sterile water for injection or 0.9% NaCl (saline)
Typical Stock Conc.	10 mg/mL (reconstitute 100 mg in 10 mL solvent)
Working Dilution	0.1 – 100 µg/mL (application-dependent)
Filter Sterilisation	0.22 µm cellulose acetate filter recommended for cell-based assays
Stability (solution)	Use within 24 hours at 2 – 8 °C; do not freeze reconstituted solution
Avoid	Calcium-containing solutions (e.g., Lactated Ringer's) — risk of precipitation

5. Recommended Research Applications

Antibacterial Testing	MIC/MBC determination; disk diffusion; broth microdilution assays
Cell Wall Biology	Studying peptidoglycan synthesis inhibition; bacterial morphology studies
Biofilm Studies	Biofilm disruption and inhibition against clinical isolates
Resistance Research	β-lactamase activity; resistance gene expression analysis
Microbiome Studies	Selective pressure assays; dysbiosis modelling
In Vitro Cell Culture	Antibiotic protection assays; intracellular infection models
NOT suitable for	In vivo (animal) or human therapeutic use

6. Handling, Safety & Storage

PPE Required	Lab coat, nitrile gloves, safety glasses; use in a well-ventilated area or fume hood
---------------------	--

Sensitisation	Known β -lactam sensitiser — avoid inhalation or skin contact; penicillin-allergic individuals should take particular care
First Aid — Skin	Wash thoroughly with soap and water
First Aid — Eyes	Rinse with copious water for ≥ 15 min; seek medical attention
First Aid — Ingestion	Do not induce vomiting; seek immediate medical advice
Disposal	Inactivate by autoclaving (121 °C, 20 min) before disposal; comply with local regulations
Storage	2 – 8 °C, in original sealed container, away from direct light and moisture
Shelf Life	24 months from date of manufacture when stored as recommended

7. Quality & Purity Specifications

Purity (HPLC)	$\geq 98\%$
Identification	IR, NMR (conforming to reference spectrum)
Loss on Drying	8.0 – 11.0% (105 °C, 4 h)
Heavy Metals	≤ 20 ppm
Residual Solvents	Within ICH Q3C limits
Microbiological	Total aerobic count ≤ 100 CFU/g
Certificate	Certificate of Analysis (COA) available upon request

8. Regulatory & Compliance Information

GHS Classification	Not classified as hazardous under GHS/WHS Regulations (Australia)
SUSMP Schedule	Schedule 4 (Prescription Only) — this product is for research use only
ADG Code	Not classified as Dangerous Goods for transport purposes
SDS Availability	Safety Data Sheet available at www.ausamics.com.au or on request
Use Restriction	For in vitro research use only by qualified personnel in licensed laboratory facilities

9. Key References

1. Neu, H.C. (1985). Cephalosporins in the treatment of serious infections. *Drug Intelligence & Clinical Pharmacy*, 19(7-8), 500–506.
2. Paladino, J.A. (1991). Ceftriaxone sodium: a third-generation cephalosporin. *Drug Intelligence & Clinical Pharmacy*, 25(10), 1075–1084.
3. Clinical and Laboratory Standards Institute (CLSI). *Performance Standards for Antimicrobial Susceptibility Testing*. Current edition.

4. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Breakpoint tables for interpretation of MICs. Current edition.

DISCLAIMER & LIABILITY NOTICE

This product is supplied for research and laboratory use only. AuSaMicS Pty Ltd (ABN 56 676 640 467) makes no representation or warranty, express or implied, regarding the fitness of this product for any particular purpose beyond its stated research application. This product is NOT approved, licensed, or intended for use in humans, animals, food processing, or clinical/diagnostic procedures. It is the sole responsibility of the end user to comply with all applicable local, state, and federal regulations governing the purchase, storage, handling, use, and disposal of this product. AuSaMicS Pty Ltd shall not be held liable for any direct, indirect, incidental, or consequential damages arising from the improper use of this product. By using this product, the purchaser agrees to these terms.