



PRODUCT DESCRIPTION

AS-1268 | LB Agar (Lennox)
www.ausamics.com.au

LB Agar (Lennox)

Reduced-Salt Bacterial Culture Medium for Molecular Microbiology

Catalogue Number: AS-1268

| Application | Medium Type | Typical pH | HS Code |
|-------------------------------|---------------------------------|--------------------|------------|
| E. coli cloning & cultivation | Non-selective solid agar medium | 7.0 ± 0.2 at 25 °C | 3821.00.00 |

Overview

LB Agar (Lennox) is a nutritionally rich, non-selective solid culture medium developed by Lennox in 1955 for the growth and maintenance of *Escherichia coli* in molecular microbiology applications. It is particularly suited to recombinant *E. coli* strains, including K-12 derivatives, which require enriched media for optimal growth and plasmid maintenance.

Compared to LB Agar (Miller), the Lennox formulation contains reduced sodium chloride at 5.0 g/L — half the NaCl of the Miller formulation (10.0 g/L) — providing lower osmolarity and improved control of osmotic conditions. This makes it well suited to salt-sensitive strains and applications where lower ionic strength supports optimal growth, improved transformation efficiency, and stable plasmid propagation.

Principle of the Medium

| | |
|----------------------------------|---|
| Tryptone (10.0 g/L) | Pancreatic digest of casein providing peptides, amino acids, nitrogen, and carbon for rapid <i>E. coli</i> growth |
| Yeast Extract (5.0 g/L) | Supplies B-vitamins, trace elements, nucleotides, coenzymes, and growth factors essential for <i>E. coli</i> K-12 |
| Sodium Chloride (5.0 g/L) | Reduced-salt Lennox formulation — lower osmolarity than Miller for salt-sensitive strains and transformation work |
| Agar (15.0 g/L) | Solidifying agent providing a firm, clear growth surface suitable for colony isolation and antibiotic selection |

Typical Composition (per litre)

| Component | CAS Number | Function | Amount |
|-----------|------------|----------|--------|
|-----------|------------|----------|--------|



PRODUCT DESCRIPTION

AS-1268 | LB Agar (Lennox)
www.ausamics.com.au

| | | | |
|--|------------|---------------------------------------|--------|
| Tryptone (Pancreatic Digest of Casein) | 73049-73-7 | Nitrogen, carbon, amino acids | 10.0 g |
| Yeast Extract | 8013-01-2 | Vitamins, growth factors, nucleotides | 5.0 g |
| Sodium Chloride | 7647-14-5 | Osmotic balance (reduced-salt Lennox) | 5.0 g |
| Agar | 9002-18-0 | Solidifying agent | 15.0 g |

Total per litre: 35.0 g | Final pH: 7.0 ± 0.2 at 25 °C

Formulation Comparison — LB Variants

| Component | Lennox (AS-1268) | Miller (AS-1267) | Bertani |
|---------------|--|-----------------------------------|---------------------|
| Tryptone | 10.0 g/L | 10.0 g/L | 10.0 g/L |
| Yeast Extract | 5.0 g/L | 5.0 g/L | 5.0 g/L |
| NaCl | 5.0 g/L | 10.0 g/L | 0 g/L |
| Agar | 15.0 g/L | 15.0 g/L | 15.0 g/L |
| Total solids | 35.0 g/L | 40.0 g/L | 30.0 g/L |
| Best for | Salt-sensitive strains, transformation | Standard E. coli — global default | Very low osmolarity |

Applications

Molecular Biology & Cloning

- Routine cultivation and maintenance of recombinant E. coli K-12 strains
- Plasmid propagation, transformation recovery, and colony isolation
- Antibiotic selection plates (ampicillin, kanamycin, chloramphenicol, tetracycline)
- Blue–white screening with X-gal and IPTG supplementation
- Salt-sensitive strain culture where lower ionic strength is required

Osmolarity-Controlled Applications

- Growth of E. coli K-12 derivatives sensitive to high NaCl concentrations
- Phage lambda infection studies — reduced NaCl is known to favour lambda propagation



PRODUCT DESCRIPTION

AS-1268 | LB Agar (Lennox)
www.ausamics.com.au

- Applications requiring precise osmotic control and reproducible transformation efficiency

Teaching & General Laboratory Use

- Standard medium for molecular biology and microbiology instruction
- General bacterial cultivation, streak plates, and competent cell preparation

Preparation Instructions

1. Suspend 35.0 g of dehydrated medium in 1 litre of demineralised or purified water.
2. Mix thoroughly and heat with agitation until fully dissolved.
3. Verify and adjust pH to 7.0 ± 0.2 if required using 1 M NaOH or 1 M HCl.
4. Sterilise by autoclaving at 121 °C for 15 minutes.
5. Cool to 45–50 °C.
6. Add filter-sterilised antibiotics or supplements aseptically as required.
7. Dispense 20–25 mL per 90 mm sterile Petri dish.
8. Allow to solidify on a level surface at room temperature before use.

Note: Do not overheat. Add heat-sensitive supplements (antibiotics, X-gal, IPTG) only after cooling to 45–50 °C.

Common Antibiotic Selection Concentrations

| Antibiotic | Working Concentration | Stock Concentration | Storage |
|-----------------|-----------------------|-----------------------|----------------------------|
| Ampicillin | 100 µg/mL | 100 mg/mL in water | -20 °C |
| Kanamycin | 50 µg/mL | 50 mg/mL in water | -20 °C |
| Chloramphenicol | 25 µg/mL | 25 mg/mL in ethanol | -20 °C |
| Tetracycline | 12.5 µg/mL | 12.5 mg/mL in ethanol | -20 °C, protect from light |
| Gentamicin | 10 µg/mL | 10 mg/mL in water | -20 °C |
| Spectinomycin | 50 µg/mL | 50 mg/mL in water | -20 °C |

Incubation & Typical Growth

Inoculate by surface spread or pour plate. Incubate aerobically at 37 °C for 12–18 hours. Colonies are visible within 12–16 hours.



PRODUCT DESCRIPTION

AS-1268 | LB Agar (Lennox)
www.ausamics.com.au

| Organism | Expected Result | Incubation |
|-----------------------------------|--|----------------|
| E. coli K-12 DH5α | Good growth — white/cream 1–3 mm colonies | 37 °C, 12–16 h |
| E. coli BL21(DE3) | Good growth — protein expression host | 37 °C, 12–16 h |
| E. coli TOP10 | Good growth — high-efficiency cloning host | 37 °C, 12–16 h |
| E. coli JM109 | Good growth — compatible with M13/phage work | 37 °C, 12–16 h |
| Salmonella typhimurium ATCC 14028 | Good growth | 37 °C, 18–24 h |

Storage & Stability

| | |
|---|--|
| Dehydrated powder | Store at 15–30 °C in a cool, dry, tightly closed container |
| Protect from | Moisture, direct light, and extreme temperatures |
| Prepared plates (no antibiotics) | 2–8 °C, inverted, up to 4 weeks |
| Prepared plates (with antibiotics) | 2–8 °C, use within 1–2 weeks |
| Shelf life (dehydrated) | As per labelled expiry date |
| Do not use if | Contamination, colour change, or cracking of agar surface observed |

Quality & Compliance

Manufactured under controlled conditions for batch-to-batch consistency. Performance tested against ATCC reference organisms prior to release. Compatible with protocols described in Lennox (1955) *Virology* 1:190, Sambrook & Russell (Molecular Cloning), and Current Protocols in Molecular Biology.

Reference

Lennox, E.S. (1955). Transduction of linked genetic characters of the host by bacteriophage P1. *Virology*, 1:190–206.

Customs & Trade Information



PRODUCT DESCRIPTION

AS-1268 | LB Agar (Lennox)
www.ausamics.com.au

HS / AHECC Code

3821.00.00

Description

Prepared culture media for the development or maintenance of microorganisms

Country of Origin

Australia

Disclaimer

This product is manufactured and supplied by AuSaMicS Pty Ltd for laboratory and research use only. It is not intended for human or veterinary consumption, therapeutic use, or in vitro diagnostic procedures without appropriate validation. AuSaMicS Pty Ltd makes no representations or warranties, express or implied, regarding the fitness of this product for any particular purpose beyond its stated intended use. Users are solely responsible for compliance with all applicable laws, regulations, and safety requirements. AuSaMicS Pty Ltd shall not be liable for any direct, indirect, incidental, or consequential damages arising from the use or misuse of this product. All information provided is believed to be accurate at the time of publication and is subject to change without notice.