



TECHNICAL DATA SHEET

AS-1364 | Thioglycollate Medium, Brewer Modified
www.ausamics.com.au

Thioglycollate Medium, Brewer Modified — Technical Data Sheet

Catalogue Number: AS-1364

Product Identification

Product Name	Thioglycollate Medium, Brewer Modified
Catalogue Number	AS-1364
HS Code	3821.00.00
Medium Type	Enriched reducing semi-solid broth for anaerobic and microaerophilic cultivation
Physical Form	Dehydrated powder
Colour (powder)	Cream to light beige with possible faint pink tint from resazurin
Colour (prepared)	Straw/amber (reduced, colourless-pink) — pink at surface indicates oxygen penetration
Pharmacopoeial equivalents	Fluid Thioglycollate Medium (USP <71>), BP 2023, EP 2.6.1
Also known as	FTM, Fluid Thioglycollate Medium, THIO Broth, Brewer Thioglycollate

Composition (per litre of prepared medium)

Ingredient	CAS Number	Function	Amount (g/L)
Pancreatic Digest of Casein	73049-73-7	Nitrogen, peptides, amino acids	15.0
Yeast Extract	8013-01-2	Vitamins, growth factors	5.0
Dextrose (D-Glucose)	50-99-7	Carbon source, mild reductant	2.5
Sodium Chloride	7647-14-5	Osmotic balance	5.0
Sodium Thioglycollate	367-51-1	Primary reducing agent	0.5
L-Cystine	56-89-3	Secondary reducing agent	0.5
Resazurin	62758-13-8	Redox indicator	0.001



TECHNICAL DATA SHEET

AS-1364 | Thioglycollate Medium, Brewer Modified
www.ausamics.com.au

Agar (semi-solid)	9002-18-0	Oxygen gradient maintenance	0.5
-------------------	-----------	-----------------------------	-----

Total: ~28.5 g/L | Final pH: 7.1 ± 0.2 at 25 °C

Physical & Chemical Specifications

Appearance (powder)	Cream to light beige homogeneous powder
Appearance (prepared)	Clear amber/straw — colourless (reduced) to faintly pink (surface only)
pH (prepared medium, 25 °C)	7.1 ± 0.2
Loss on Drying (moisture)	≤ 5.0%
Solubility	Completely soluble with gentle heating
Agar concentration	0.5 g/L (semi-solid — not solid agar)
Resazurin	0.001 g/L — colourless (reduced/anaerobic) ↔ pink (oxidised/aerobic)
Reducing capacity	Dual system: sodium thioglycollate + L-cystine
Pharmacopoeial grade	Meets USP <71>, BP 2023, EP 2.6.1 Fluid Thioglycollate Medium requirements

Performance Specifications

Organism	ATCC No.	Incubation	Expected Result	Pharmacopoeial Ref.
Clostridium sporogenes	11437	30–35 °C, 5 days (anaerobic)	Growth in lower zone — obligate anaerobe	USP <71> / EP 2.6.1
Bacteroides vulgatus	8482	35–37 °C, 48–72 h (anaerobic)	Growth in lower zone	USP <71>
Staphylococcus aureus	6538	30–35 °C, 3 days (aerobic)	Growth throughout	USP <71> / EP 2.6.1
Pseudomonas aeruginosa	9027	30–35 °C, 3 days	Growth at surface	USP <71> / EP 2.6.1
Candida albicans	10231	20–25 °C, 5 days	Growth throughout	USP <71>



TECHNICAL DATA SHEET

AS-1364 | Thioglycollate Medium, Brewer Modified
www.ausamics.com.au

Bacillus subtilis	6633	30–35 °C, 5 days	Growth at surface	USP <71> / EP 2.6.1
-------------------	------	------------------	-------------------	---------------------

Preparation Protocol

1. Dissolve 28.5 g per litre in purified water with gentle heating and agitation.
2. Autoclave at 121 °C for 15 minutes.
3. Cool to room temperature.
4. Before use: heat in boiling water bath or steamer for 10–15 minutes. Cool rapidly without shaking.
5. Discard any tube where >1/3 of depth is pink (indicating oxidation).

For pharmacopoeial sterility testing: follow the specific preparation, validation, and incubation requirements of the applicable pharmacopoeia (USP, BP, or EP).

Key Limitations

- Presumptive medium for anaerobic detection — confirmatory subculture required for organism identification
- Resazurin pink layer must not exceed 1/3 of tube depth before inoculation
- Not suitable as the sole medium for definitive anaerobe identification
- For USP/BP/EP sterility testing: medium must be prepared and validated per pharmacopoeial requirements
- Do not overheat — excessive autoclaving degrades reducing agents and nutrients

Storage & Stability

Dehydrated powder	15–30 °C, tightly sealed, dry
Prepared medium	2–8 °C, protected from light, minimise oxygen exposure
Shelf life	As per labelled expiry
Before use	Heat to expel dissolved O ₂ ; discard if >1/3 pink

Regulatory & Quality Standards

Pharmacopoeial compliance	USP <71> — Sterility Tests; BP 2023; EP 2.6.1
Country of Manufacture	Australia



TECHNICAL DATA SHEET

AS-1364 | Thioglycollate Medium, Brewer Modified
www.ausamics.com.au

Performance testing

Per USP <71> / EP 2.6.1 reference organisms

Disclaimer

This product is manufactured and supplied by AuSaMics Pty Ltd for laboratory and research use only. It is not intended for human or veterinary consumption, therapeutic use, or in vitro diagnostic procedures without appropriate validation. AuSaMics Pty Ltd makes no representations or warranties, express or implied, regarding the fitness of this product for any particular purpose beyond its stated intended use. Users are solely responsible for compliance with all applicable laws, regulations, and safety requirements. AuSaMics Pty Ltd shall not be liable for any direct, indirect, incidental, or consequential damages arising from the use or misuse of this product. All information provided is believed to be accurate at the time of publication and is subject to change without notice.