



PRODUCT DESCRIPTION

AS-1395 | 2× YT Agar
www.ausamics.com.au

2× YT Agar

Double-Strength Yeast Tryptone Agar — Rich Medium for High-Density *E. coli* Growth
Catalogue Number: AS-1395

Application	Medium Type	Typical pH	HS Code
High-density <i>E. coli</i> & recombinant	Non-selective solid agar medium	7.0 ± 0.2 at 25 °C	3821.00.00

Overview

2× YT Agar (Double-Strength Yeast Tryptone Agar) is a nutrient-rich, non-selective solid culture medium designed to support rapid and high-density growth of bacteria, particularly *Escherichia coli* used in molecular biology and recombinant DNA workflows.

Compared to standard LB or YT media, 2× YT Agar contains higher concentrations of yeast extract (10.0 g/L) and tryptone (16.0 g/L), providing an increased supply of peptides, amino acids, vitamins, and growth factors. This enhanced formulation results in larger colonies, faster growth, and improved plasmid yield in cloning, expression, and phage-related applications.

Comparison — 2× YT vs LB vs SOB

Component	2× YT (AS-1395)	LB Miller (AS-1267)	SOB
Tryptone	16.0 g/L	10.0 g/L	20.0 g/L
Yeast Extract	10.0 g/L	5.0 g/L	5.0 g/L
NaCl	5.0 g/L	10.0 g/L	0.5 g/L
Agar	15.0 g/L	15.0 g/L	15.0 g/L
Total solids	46.0 g/L	40.0 g/L	40.5 g/L
Best for	High-density, fast growth	Standard routine	High-efficiency transformation

Principle of the Medium



PRODUCT DESCRIPTION

AS-1395 | 2× YT Agar
www.ausamics.com.au

Tryptone (16.0 g/L)	Pancreatic digest of casein — high-concentration source of peptides, amino acids, nitrogen, and carbon for rapid E. coli growth
Yeast Extract (10.0 g/L)	Double-strength supply of B-vitamins, trace elements, nucleotides, coenzymes, and growth factors to sustain high cell density
Sodium Chloride (5.0 g/L)	Maintains osmotic balance compatible with most laboratory E. coli strains
Agar (15.0 g/L)	Solidifying agent providing a firm, clear growth surface for colony isolation and antibiotic selection

Typical Composition (per litre)

Component	CAS Number	Function	Amount
Tryptone (Pancreatic Digest of Casein)	73049-73-7	High-level nitrogen, carbon, amino acids	16.0 g
Yeast Extract	8013-01-2	Double-strength vitamins, growth factors	10.0 g
Sodium Chloride	7647-14-5	Osmotic balance	5.0 g
Agar	9002-18-0	Solidifying agent	15.0 g

Total per litre: 46.0 g | Final pH: 7.0 ± 0.2 at 25 °C

Key Features & Benefits

- Double-strength tryptone and yeast extract for enhanced bacterial growth and biomass yield
- Supports larger colonies and faster growth than standard LB medium
- Ideal for high-density E. coli culture prior to protein induction (IPTG)
- Improved plasmid yield per colony due to enriched nutrient environment
- Suitable for antibiotic-supplemented selection plates
- Compatible with all standard E. coli K-12 and B strains
- Widely used in phage display, M13 phage propagation, and panning workflows

Applications

Molecular Biology & Genetics

- High-density cultivation of E. coli for cloning and plasmid propagation



PRODUCT DESCRIPTION

AS-1395 | 2× YT Agar
www.ausamics.com.au

- Growth of recombinant and expression strains (BL21, Rosetta, ArcticExpress)
- Pre-induction culture growth prior to IPTG addition for protein expression
- Antibiotic selection plates — ampicillin, kanamycin, chloramphenicol
- Blue–white screening with X-gal and IPTG supplementation

Phage & Protein Expression Work

- Bacteriophage propagation — M13, lambda, T4 phage
- High-density bacterial lawns for plaque assays
- Phage display library screening (panning on agar)
- Protein expression studies requiring high biomass

Research & Teaching Laboratories

- Routine bacterial culture where fast growth is required
- Demonstration of nutrient-rich vs minimal media growth rates
- Microbial physiology and growth kinetics studies

Preparation Instructions

1. Suspend 46.0 g of dehydrated medium in 1 litre of purified or demineralised water.
2. Heat with agitation until completely dissolved.
3. Verify and adjust pH to 7.0 ± 0.2 if required using 1 M NaOH or 1 M HCl.
4. Sterilise by autoclaving at 121 °C for 15 minutes.
5. Cool to 50–55 °C.
6. Add filter-sterilised antibiotics or supplements aseptically if required.
7. Mix gently and pour 20–25 mL per 90 mm sterile Petri dish.
8. Allow to solidify on a level surface before use.

Note: Do not overheat. The slightly higher solidification temperature (50–55 °C vs 45–50 °C for LB) is recommended to prevent premature solidification. Add heat-sensitive supplements only after cooling.

Common Antibiotic Selection Concentrations

Antibiotic	Working Concentration	Stock Concentration	Storage
Ampicillin	100 µg/mL	100 mg/mL in water	–20 °C
Kanamycin	50 µg/mL	50 mg/mL in water	–20 °C



PRODUCT DESCRIPTION

AS-1395 | 2× YT Agar
www.ausamics.com.au

Chloramphenicol	25 µg/mL	25 mg/mL in ethanol	-20 °C
Tetracycline	12.5 µg/mL	12.5 mg/mL in ethanol	-20 °C, protect from light
Gentamicin	10 µg/mL	10 mg/mL in water	-20 °C

Incubation & Typical Growth

Inoculate by surface spread or pour plate. Incubate aerobically at 37 °C for 12–18 hours. Colonies on 2× YT Agar are typically larger and more vigorous than on LB medium due to the enriched nutrient content.

Organism	Expected Result	Incubation
E. coli K-12 DH5α	Good growth — large cream colonies, >3 mm	37 °C, 12–16 h
E. coli BL21(DE3)	Vigorous growth — suitable for expression work	37 °C, 12–16 h
E. coli TOP10	Good growth — high-efficiency cloning	37 °C, 12–16 h
E. coli XL1-Blue	Good growth — M13 phage compatible	37 °C, 12–16 h
Salmonella typhimurium ATCC 14028	Good growth	37 °C, 18–24 h

Storage & Stability

Dehydrated powder	Store at 15–30 °C in a cool, dry, tightly closed container
Protect from	Moisture, direct light, and extreme temperatures
Prepared plates (no antibiotics)	2–8 °C, inverted, up to 4 weeks
Prepared plates (with antibiotics)	2–8 °C, use within 1–2 weeks
Shelf life (dehydrated)	As per labelled expiry date
Do not use if	Contamination, colour change, or cracking of agar surface observed

Quality & Compliance



PRODUCT DESCRIPTION

AS-1395 | 2× YT Agar
www.ausamics.com.au

Manufactured under controlled conditions for batch-to-batch consistency. Performance tested against ATCC reference organisms prior to release. Compatible with protocols described in Sambrook & Russell (Molecular Cloning), Current Protocols in Molecular Biology, and phage display methodologies (Barbas et al., Phage Display: A Laboratory Manual, CSHL Press).

Customs & Trade Information

HS / AHECC Code	3821.00.00
Description	Prepared culture media for the development or maintenance of microorganisms
Country of Origin	Australia

Disclaimer

This product is manufactured and supplied by AuSaMicS Pty Ltd for laboratory and research use only. It is not intended for human or veterinary consumption, therapeutic use, or in vitro diagnostic procedures without appropriate validation. AuSaMicS Pty Ltd makes no representations or warranties, express or implied, regarding the fitness of this product for any particular purpose beyond its stated intended use. Users are solely responsible for compliance with all applicable laws, regulations, and safety requirements. AuSaMicS Pty Ltd shall not be liable for any direct, indirect, incidental, or consequential damages arising from the use or misuse of this product. All information provided is believed to be accurate at the time of publication and is subject to change without notice.