

Ascospore Agar

Catalogue No. AS-1436 · Dehydrated Culture Medium

Section 1 — Product Identification

Product Name	Ascospore Agar
Catalogue Number	AS-1436
Product Format	Dehydrated powder
Intended Use	Selective detection and isolation of ascosporegenous yeasts from food, beverage, environmental, and industrial samples
Manufacturer	AuSaMics Pty Ltd
Address	31 Longview CT, Thomastown, VIC 3074, Australia
ABN	56 676 640 467
Website	www.ausamics.com.au
Email	support@ausamics.com.au
Phone	+61 412 520 598

Section 2 — Mode of Action

Ascospore Agar is a selective medium formulated to promote ascosporegenesis in yeast species. Potassium acetate serves as the primary carbon source and acts as a sporulation stimulus by creating a nutritionally stressful environment that triggers yeast cells to form ascospores. Yeast extract and dextrose provide minimal nutrition to maintain cell viability during incubation. The high agar concentration (30 g/L) produces a firm medium that supports colony morphology examination. The selective pressure imposed by limited nutrients inhibits non-sporulating yeasts and most bacteria, enabling the detection of Zygosaccharomyces, Saccharomyces, and related ascosporegenous genera from processed food and beverage samples.

Section 3 — Composition

Ingredient	g / L
Yeast Extract	2.5
Dextrose (Glucose)	1.0
Potassium Acetate	10.0
Agar	30.0
Total (per litre)	43.5 g

Final pH at 25°C: 6.4 ± 0.2

Section 4 — Preparation

1. Suspend 43.5 g of the dehydrated powder in 1 litre of distilled or deionised water.
2. Mix thoroughly and heat to boiling with agitation until completely dissolved.
3. Sterilize by autoclaving at 121°C for 15 minutes.
4. Allow to cool to 45–50°C before pouring into Petri dishes.
5. Allow plates to solidify and dry before use.

Note: Do not overheat — excessive heating may reduce gel strength and alter pH.

Section 5 — Quality Control

Organism (ATCC)	Inoculum (CFU)	Incubation	Expected Result
<i>Saccharomyces cerevisiae</i> (9763)	≤100	25±2°C / 3–5 days	Good growth · ascospore formation
<i>Zygosaccharomyces rouxii</i> (2623)	≤100	25±2°C / 3–5 days	Good growth · ascospore formation
<i>Escherichia coli</i> (25922)	≤100	25±2°C / 3–5 days	Inhibited / No growth
<i>Staphylococcus aureus</i> (25923)	≤100	25±2°C / 3–5 days	Inhibited / No growth

Dehydrated Appearance: Cream-white to light beige, free-flowing homogeneous powder.

Prepared Appearance: Light amber, clear to slightly opalescent gel.

Section 6 — Storage

Dehydrated Powder Store at 15–25°C in a tightly sealed container, away from moisture and direct sunlight

Prepared Medium Store at 2–8°C; use within 4 weeks of preparation

Shelf Life Refer to expiry date on label

Section 7 — Applicable Standards

Standard	Relevance
ISO 11133:2014	Culture media performance testing — QC criteria
ISO 21527-2:2008	Yeasts and moulds in food with low water activity (osmotolerant)
ISO 7954:1987	Enumeration of yeasts and moulds (general method)

Section 8 — Literature / References

#	Reference	Relevance
1	Lodder J, Kreger-van Rij NJW. The Yeasts: A Taxonomic Study. Amsterdam: North-Holland; 1952.	Classical taxonomy of ascosporegenous yeasts
2	International Organization for Standardization. ISO 21527-2:2008. Microbiology of food and animal feeding stuffs — Horizontal method for the enumeration of yeasts and moulds — Part 2: Colony count technique in products with water activity less than or equal to 0.95. Geneva: ISO; 2008.	Primary standard for osmotolerant yeast enumeration
3	International Organization for Standardization. ISO 11133:2014. Microbiology of food, animal feed and water — Preparation, production, storage and performance testing of culture media. Geneva: ISO; 2014.	QC performance criteria
4	Beuchat LR. Thermal inactivation of yeasts in fruit juices supplemented with food-grade additives. J Food Sci. 1981;46(4):1137–1139.	Heat-resistant yeast detection in food

Disclaimer: The information provided in this document is accurate to the best of AuSaMicS Pty Ltd's knowledge at the time of publication. It is intended solely as a guide for the safe and effective use of this product. AuSaMicS Pty Ltd accepts no liability for any loss, damage, or injury arising from the use of this product outside its intended application. Users are responsible for ensuring suitability for their specific application. This document does not constitute a warranty of any kind. All regulatory obligations remain the responsibility of the end user.