



## Technical Data Sheet

TDS-AS-1151  
www.ausamics.com.au

### Product Identification

Product Name	Brilliant Green Bile Agar (BGBA)
Catalog No.	AS-1151
Lot No.	BGB260301
Mfg. Date	March 2026
Retest Date	March 2028
Grade	Microbiological Grade — Highly Selective Salmonella Medium
Formulation	Based on Kristensen/Kauffmann Brilliant Green Bile formulation

### Physical & Chemical Properties

Appearance (powder)	Homogeneous; light beige to pale yellow-green powder
Appearance (prepared)	Green-amber gel; medium turns brownish-amber on cooling
Odour	Faint, characteristic (peptone/bile)
pH (prepared, 25 °C)	7.2 ± 0.2
Agar gel / melt	Gels ~32–34 °C / Melts ~84–86 °C
Loss on Drying	≤ 7.0% w/w (105 °C / 2 h)
Suspension	55.1 g / L — dissolves completely on boiling (NO AUTOCLAVE)

### Composition per Litre

Component	Function	Amount (g/L)
Peptone	Nitrogen, amino acids — supports Salmonella colony size	10.0
Meat Extract	Growth factors, vitamins, minerals	5.0



## Technical Data Sheet

TDS-AS-1151  
www.ausamics.com.au

Component	Function	Amount (g/L)
Ox Bile (dried) — 2%	Primary selective agent — near-complete suppression of coliforms and Proteus	20.0
Sodium Chloride	Osmotic balance	5.0
Brilliant Green	Inhibits Gram-positive and most coliform Gram-negatives; synergises with bile	0.0047
Phenol Red	pH indicator — Salmonella: pink-red; no fermentable sugars so no colour masking	0.09
Agar	Solidifying agent	15.0
TOTAL		55.1

**CRITICAL: DO NOT AUTOCLAVE — Boiling only. Autoclaving destroys brilliant green and bile salt selectivity.**

### Preparation Protocol

Step 1 — Suspension	Suspend 55.1 g of dehydrated BGBA in 1 L of purified water.
Step 2 — Dissolution	Heat with frequent agitation. Bring to a full boil for 1 minute until completely dissolved.
Step 3 — NO AUTOCLAVE	DO NOT autoclave. Autoclave destroys brilliant green and bile salts, eliminating selectivity and differential properties.
Step 4 — Cooling	Cool to 45–50 °C with gentle mixing.
Step 5 — Pouring	Pour 15–20 mL per 90 mm Petri dish. Allow to solidify on a level surface.
Step 6 — Drying	Dry plates at 37 °C for 15–20 min. Store at 2–8 °C, protected from light. Use within 7 days.
Inoculation	Surface-plate 0.1–0.5 mL enrichment broth. Incubate 35–37 °C for 18–24 h.

### Incubation & Colony Reading

Temperature	35–37 °C
-------------	----------



## Technical Data Sheet

TDS-AS-1151  
www.ausamics.com.au

Duration	18–24 h primary; up to 48 h for slow growers
Salmonella	Pink-red to red colonies; 2–5 mm; smooth, convex, well-defined
Coliforms	Completely or almost completely inhibited
Proteus	Completely inhibited
Confirmation	All suspect colonies must be confirmed by TSI, API 20E, VITEK, PCR, and/or serotyping
Note	Not suitable for <i>S. Typhi</i> and <i>S. Paratyphi</i>

### Quality Control Specifications

Parameter	Specification	Method
Appearance (powder)	Homogeneous; light beige to pale yellow-green; no lumps	Visual inspection
pH (prepared, 25 °C)	7.0 – 7.4	Potentiometry (ISO 10523)
Loss on Drying	≤ 7.0% w/w	105 °C / 2 h (gravimetric)
Growth — <i>S. typhimurium</i> ATCC 14028	Good growth; pink-red 2–5 mm; ≥ 50% recovery	ISO 11133 growth test
Growth — <i>S. enteritidis</i> ATCC 13076	Good growth; pink-red colonies; ≥ 50% recovery	ISO 11133 growth test
Inhibition — <i>E. coli</i> ATCC 25922	Complete / ≥ 99% inhibition	ISO 11133 selectivity
Inhibition — <i>Proteus mirabilis</i> ATCC 12453	Complete / ≥ 99% inhibition	ISO 11133 selectivity
Sterility (prepared medium)	No growth at 14 days (30–35 °C)	Incubation sterility check

### Limitations

Typhoidal Salmonella	NOT suitable for <i>S. Typhi</i> and <i>S. Paratyphi</i> — use Bismuth Sulfite or SS Agar
Presumptive Only	BGBA is a presumptive medium; all suspect colonies require biochemical and serological confirmation



## Technical Data Sheet

TDS-AS-1151  
www.ausamics.com.au

No H <sub>2</sub> S indicator	Does not differentiate H <sub>2</sub> S-producing Salmonella; use alongside XLD or BS Agar
Autoclave Warning	Autoclaving destroys selectivity — product must be prepared by boiling only

### Technical References

Original Reference	Kristensen M., Lester V. & Jorgens A. (1925). <i>Über die Verwendung von tetrathionat-nährboden und von brillant-grün bei der Typhusdiagnose.</i>
Modification	Kauffmann F. & Kristensen M. — Brilliant Green Bile Agar, modified formulation
ISO Standard	ISO 6579-1:2017 — Horizontal method for detection of Salmonella spp.
FDA-BAM	FDA Bacteriological Analytical Manual Chapter 5 — Salmonella
Historical Use	Referenced in German and Nordic national food microbiology standards

**Disclaimer:** *AuSaMicS Pty Ltd (ABN 56 676 640 467) warrants that this product meets the specifications stated herein at time of manufacture and release. This product is for laboratory use only. Not intended for human, veterinary, or food consumption. AuSaMicS Pty Ltd shall not be held liable for any direct, indirect, incidental, or consequential damages resulting from use outside its intended application. By using this product the purchaser agrees to these terms.*